INSTRUCTIONS FOR USE

INTENDED USE

CD-Chex TdT™ Plus is intended to be used as a quality control material for evaluating intracellular and surface antigens, including TdT (Terminal deoxynucleotidyl transferase), CD1a, CD34 and Cytoplasmic CD3, with monoclonal antibody binding by flow cytometry. When these cells are stained with fluorescent antibodies and analyzed by flow cytometry, they provide a reference value for abnormal cells found in certain types of hematopoietic neoplasms. CD-Chex TdT Plus is designed for use on BD Biosciences and Beckman Coulter® flow cytometry systems. **This product and the markers provided on the assay have not been cleared by the U.S. Food and Drug Administration for In Vitro Diagnostic use. This product and the values provided are For Research Use Only. Not for use in diagnostic procedures.**

SUMMARY AND PRINCIPLES

Immunophenotyping by flow cytometry is a complex process. CD-Chex TdT Plus is designed to represent abnormal peripheral blood leukocytes consistent with a hematolymphoid neoplastic patient sample<sup>1,2,3</sup>. CD-Chex TdT Plus mimics whole blood samples by possessing surface antigens and intracellular antigens detectable with fluorescent monoclonal antibodies by flow cytometry. Stable antigens on CD-Chex TdT Plus include surface CD1a, CD2, CD4, CD5, CD7, CD8, CD34, intracellular CD3 and nuclear TdT. CD-Chex TdT Plus is a positive procedural assayed control used to monitor reagent staining, erythrocyte lysis, sample preparation, and instrument performance to provide consistent and reliable quality control measurements<sup>4</sup>.

REAGENTS

CD-Chex TdT Plus contains stabilized human blood and cells of human origin in a preservative medium.

PRECAUTIONS

1. CD-Chex TdT Plus is For Research Use Only. Not for use in diagnostic procedures.
2. All human source material used to manufacture this product was non-reactive for antigens to Hepatitis B (HBsAg), negative by tests for antibodies to HIV (HIV-1/HIV-2) and Hepatitis C (HCV), non-reactive for HIV-1 RNA and HCV RNA by licensed NAT, and non-reactive to Serological Test for Syphilis (STS) using techniques specified by the U.S. Food and Drug Administration. Because no known test method can assure complete absence of human pathogens, this product should be handled with appropriate precautions.
3. This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
4. This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product as supplied invalidates the intended use of the product.
5. CD-Chex TdT Plus should not to be used as calibrator.

STORAGE AND STABILITY

CD-Chex TdT Plus is stable through the expiration date when stored at 2 °C to 10 °C. After opening, the product is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C. DO NOT FREEZE.

INDICATIONS OF PRODUCT DETERIORATION

Inability to obtain expected values may indicate product deterioration. If the recovered values are not within the expected ranges:

1. Review control product package insert and the operating procedure of the instrument.
2. Check expiration date of the product on the vial. Discard outdated products.
3. Assay an unopened vial of the CD-Chex TdT Plus. If the values are still outside the expected range, contact Streck Technical Services at 800-843-0912 or [technicalservices@streck.com](mailto:technicalservices@streck.com).

INSTRUCTIONS FOR USE

1. Follow instrument manufacturer's instructions for instrument compensation and sample analysis.
2. Remove a vial of the control from refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
3. To mix:
  - a. Hold horizontally between palms of the hands and roll the vial back and forth for 20 to 30 seconds.
  - b. Mix by gentle inversion at least 8-10 times until all cells are thoroughly suspended.
  - c. Gently invert the vial 5 times immediately before sampling.

*Note: Vials stored for an extended period of time may require extra mixing.*

4. Aliquot appropriate volume of control reagent into each tube and return to refrigeration to ensure maximum open vial stability.
5. Add recommended monoclonal antibodies according to manufacturer's instructions to each tube and mix gently.

*Note: A negative staining control is recommended due to heterogeneous/dim expression of selected assayed parameters.*

6. Incubate according to antibody manufacturer's instructions.
7. Add recommended amount of RBC lysing agent according to manufacturer's instructions.
8. Analyze by flow cytometry in the same manner as patient samples.

INTRACELLULAR MARKER INSTRUCTIONS FOR USE

Follow steps 1-3 above. CD-Chex TdT Plus is stabilized cellular material. Therefore, the fixation step (Reagent 1 or Reagent A) used prior to permeabilization in commercially available intracellular staining kits is not required. Use of the fixation reagent will result in sub-optimal recovery.

LIMITATIONS

CD-Chex TdT Plus is not intended to be used in ISHAGE protocols for CD34 enumeration. Staining characteristics are not equivalent to the phenotypic properties of human hematopoietic progenitor cells<sup>5,6</sup>.

EXPECTED RESULTS

The mean values provided are derived from replicate analyses on properly compensated flow cytometers. The assay values are obtained using reagents recommended by each instrument manufacturer. The expected ranges listed represent estimates of variation due to different laboratories' reagents, instrument performance, maintenance, and operator technique.

It is recommended that an individual laboratory establish its own control means and ranges for each parameter and sample analyzed that reflect the laboratory's specific conditions and protocols. These established control means should fall within the published expected ranges. Data collected from interlaboratory quality control programs can be used as a cumulative approach when calculating ranges.

GATING

The most common gating strategy used in neoplastic cell assessment is to gate on the abnormal cells and then determine the CD marker percent positivity using a negative staining control<sup>1-4</sup>. Abnormal cells are generally located on a FSC/SSC plot or a CD45/SSC plot; although other gating strategies can be employed<sup>1-4</sup>. Debris should be excluded from the cell gate to obtain percent recovery values within assay range.


FIGURE 1. Gating methods used to identify abnormal cells.

In CD-Chex TdT Plus, abnormal cells can be identified using standard gating strategies: FSC/SSC (A. Rectangle Gate) or CD45/SSC (B. Oval Gate). Optimal gates will remove debris for evaluation of percent recovery.

REFERENCES

1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press, Northfield, IL USA, 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematolymphoid cells. Approved Guideline - Second Edition, Wayne, PA USA, 2009.
5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition, Wayne, PA USA, 2007.
6. Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

QUALITY CONTROL PROGRAM

Streck offers STATS®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the STATS Department at 800-898-9563 or [statsdata@streck.com](mailto:statsdata@streck.com). Additional information can be found at [www.streck.com](http://www.streck.com).

ORDERING INFORMATION

Please call our Customer Service Department toll free 800-228-6090 for assistance. Additional information can be found online at [www.streck.com](http://www.streck.com).

GLOSSARY OF HARMONIZED SYMBOLS

Batch Code	Biological Risk	Catalog Number	Use By	Manufacturer	Consult Instructions For Use	Temperature Limitation
Glossary of symbols may contain symbols not used in the labeling of this product.						

See [www.streck.com/patents](http://www.streck.com/patents) for patents that may be applicable to this product.

The brand and product names of the instruments are trademarks of their respective holders.

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### RÉFÉRENCES

1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press, Northfield, IL USA, 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematoïmaphoid cells. Approved Guideline - Second Edition, Wayne, PA USA, 2009.
5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition, Wayne, PA USA, 2007.
6. Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

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### QUELLENANGABEN

- Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
- Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press, Northfield, IL USA, 2012.
- Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
- Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematolymphoid cells. Approved Guideline - Second Edition, Wayne, PA USA, 2009.
- Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition, Wayne, PA USA, 2007.
- Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

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### BIBLIOGRAFIA

1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press, Northfield, IL USA, 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hemolymphoid cells. Approved Guideline - Second Edition, Wayne, PA USA, 2009.
5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition, Wayne, PA USA, 2007.
6. Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherapy 1996; 5: 213-226.

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1. Gorczyca, W. Flow cytometry in neoplastic hematology. Morphologic-Immunophenotypic Correlation - Second Edition, Informa Healthcare, New York, NY USA, 2010.
2. Cheria, S., Wood, B. Flow cytometry in hematopoietic neoplasms; A case-based approach. CAP Press, Northfield, IL USA, 2012.
3. Nguyen, D., Diamond, L.W., Braylan, R.C. Flow cytometry in hematopathology, a visual approach to data analysis and interpretation. Second Edition, Humana Press, Totowa, NJ USA, 2007.
4. Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematomphoid cells. Approved Guideline - Second Edition, Wayne, PA USA, 2009.
5. Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition. Wayne, PA USA, 2007.
6. Sutherland, D.R., Anderson L., Keeney, M., Nayar, R., Chin-Yee, I. The ISHAGE guidelines for CD34+ cell determination by flow cytometry. J Hematotherap 1996; 5: 213-226.

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